

## Product Datasheet

### EBV-VCA gp125, Purified Ag (orb358022)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | EBV-VCA gp125, Purified Ag                              |
| <b>Tested Applications</b> | CLIA, EIA, LF, WB                                       |
| <b>Target</b>              | Epstein-Barr Virus (EBV-VCA) Viral Capsid Antigen gp125 |
| <b>Preservatives</b>       | No Preservative                                         |
| <b>Note</b>                | For research use only                                   |
| <b>Source</b>              | P3H3 Cells                                              |
| <b>Expiration Date</b>     | 6 months from date of receipt.                          |