

Product Datasheet

C1orf87 antibody (orb35809)

Description	C1orf87 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein C1orf87 protein (1-138AA)
Target	C1orf87
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N0U7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.