


Product Datasheet

Human CD206 protein (orb358245)

Description	The MMR (macrophage mannose receptor) is also called MR due to its presence on cells other than macr...
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	RTSLCFKLYAKGKHEKKTWFESRDFCRALGGDLASINNKEEQTIWRLITASGSY HKLFWLGLTYGSPSEGFTWSDGSPVSYENWAYGEPNNYQNVEYCGELKGDPT MSWNDINCEHLNNWICQIQKGQTPKPEPTPAPQDNPPVTEDGWVIYKDYQYYFS KEKETMDNARAFCKRNFGDLVSIQSESEKKFLWKYVNRNDAQSAYFIGLLISLDK KFAWMDGSKVDYVSWATGEPNFADEDENCVTMYSNSGFWNDINCGYPNAFIC QRHNSSINATTVMPTMPSVPSGCKEGWNFYSNKCFKIFGFMEERKNWQEAR ACIGFGGNLVSQNEKEQAFLTYHMKDSTFSAWTGLNDVNSEHTFLWTDGRGV HYTNWVGKGYPGGRRSSLSYEDADCVVIIGGASNEAGKWMDDTCDSKRGYICQT RSDPSLTNPPATIQTDFVKYKSSYSLMRQKFQWHEAETYCKLHNSLIASILDPY SNAFAWLQMETSNERVWIALNSNLTDNQYTWTDKWRVRYTNWAADEPKLKSA CVYLDLDGYWKTACNESFYFLCK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	66.4 kDa
Uniprot ID	P22897
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAG...