




Product Datasheet

BLOC1S2 antibody (orb35833)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	BLOC1S2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15961
Uniprot ID	Q6QNY1
NCBI	NP_776170.2 , NP_001001342.1 , NP_001269367.1 , NP_001269366.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of BLOC1S2 antibody...



Flow cytometric analysis of BLOC1S2 anti...