

Product Datasheet

Mouse Hspa1b protein (orb358342)

Description	Recombinant mouse Hspa1b protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	AKNTAIGIDLGGTYSVGVFQHGKVEIIANDQGNRTTPSYAFTDTERLIGDAAKN QVALNPQNTVFDAKRLIGRKFGDAVVQSDMKHWPQVVDGDKPKVQVNYKG ESRSFFPEEISSMVLTKMKEIAEAYLGHPVTNAVITVPAYFNDSQRQATKDAGVIA GLNVLRINIPTAAAIAYGLDRTGKGERNVLIFDLGGGTFDVSILTIDDGIFEVKAT AGDTHLGGEDFDNRLVSHFVEEFKRKHKKDISQNKRAVRRRLRTACERAKRTLSS STQASLEIDSLFEGIDFYTSITRARFEELCSDLFRGTLEPVEKALRDAKMDKAQIHD LVLVGGSTRIPKVQKLLQDFFNGRDLNKSINPDEAVAYGAAVQAAILMGDKSEN VQDLLLLDVAPLSLGLETAGGVMTALIKRNSTIPTKTQTFTTYSNQPGLVIQVY EGERAMTRDNNLLGRFELSGIPPAPRGVPQIEVTFDIDANGILNVTATDKSTGKA NKITITNDKGRLSKEEIERMVQEAERYKAEDEVQRDRVAAKNALESYAFNMKSAV EDEGLKGLSEADKKKVLDKCQEVISWLDSENTLADKEEFVHKREELERVCSPIIS GLYQGAGAPGAGGFQAQAPPKGASGSGPTIEEVD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	86 kDa
Uniprot ID	P17879
Expiration Date	6 months from date of receipt.