

Product Datasheet

Bacterial aprE protein (orb358396)

Description	Recombinant bacterial aprE protein
Endotoxins	Not test.
Preservatives	Tris/PBS-based buffer, 5%-50% glycerol
Form/Appearance	Liquid
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	AQSVPYGISQIKAPALHSQGYTGSNVKVAVIDSGIDSSHPDLNVRGGASFPSET NPYQDGSSSHGTHVAGTIAALNNSIGVLGVAPSASLYAVKVLDTGSGQYSWIING IEWAISNNMDVINMSLGGPTGSTALKTVVDKAVSSGIVVAAAAGNEGSSGSTST VGYPAYKPSTIAVGAVNSSNQRAFSSAGSELDVMAPGVSIQSTLPGGTYGAYN GTSMATPHVAGAAALILSKHPTWTNAQVRDRLESTATYLGNSFYFGKGLINVQA AAQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	43.7 kDa
Uniprot ID	P04189
Expiration Date	6 months from date of receipt.