

Product Datasheet

C17orf78 antibody (Biotin) (orb35848)

Description	C17orf78 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein C17orf78 protein (1-186AA)
Target	C17orf78
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N4C9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.