

Product Datasheet

Human KKLC1 protein (orb358909)

Description	Recombinant human KKLC1 protein
Endotoxins	Less than 1.0 EU/ug as determined by LAL method.
Preservatives	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged(This tag can be tested only under denaturing conditions)
Note	For research use only
Application notes	Full length of HIS-tag and expression region is 20.63kDa
Protein Sequence	MN FYLL LASSIL CALIVFWKYRRFQRNTGEMSSNSTALALVRPSSSGLINSNTDN NLAVYDL SRDILNNFPHSIARQKRILVNLSMVENKLVELEHTLLSKGFRGASPHRK ST
Purity	Testing in progress.
Source	Mammalian cell
MW	14 kDa
Uniprot ID	Q5H943
Expiration Date	6 months from date of receipt.