

Product Datasheet

Rat Lum protein (orb358934)

Description	Recombinant rat Lum protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full length of HIS-tag protein
Protein Sequence	QYYDYDAPLFMYGELSPNCAPECNCPHSYPTAMYCDDLKLSVPMVPPGIKYLYL RNNQIDHIDEKAFENVTDLQWLILDHNLENSKIKGVFSKLKQLKKLHINYNNLT ESVGPLPKSLQDLQLANNKISKLGSDGLVNLTFIYLQHNQLKEEAVSASLKGLKS LEYLDLSFNQMSKLPAGLPTSLLTLYLDNNKITNIPDEYFNRFTGLQYLRLSHNELA DSGVPGNFSNLSLLELDLSYNKLKSIPTVNENLENYYLEVNKLEKFDVKSFKILG PLSYSKIKHLRLDGNPLTQSSLPDMECLRVANEITVN
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	40.5 kDa
Uniprot ID	P51886
Expiration Date	6 months from date of receipt.