




Product Datasheet

REG3G antibody (orb35898)

Description	Rabbit polyclonal antibody to REG3G first isolated from a regenerating islet cDNA library, encodes a...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	REG3G
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19330
Uniprot ID	Q6UW15
NCBI	NP_001008388.1 , NP_940850.1 , NP_001256969.1
Dilution Range	IHC-P: 1:50-100, WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of REG3G binding...



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections.