



## Product Datasheet

### GCDH antibody (orb35911)


|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | GCDH antibody                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                         |
| <b>Reactivity</b>          | Human, Mouse                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | ELISA, IHC, IP, WB                                                             |
| <b>Immunogen</b>           | Recombinant Human Glutaryl-CoA dehydrogenase, mitochondrial protein (45-300AA) |
| <b>Target</b>              | GCDH                                                                           |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                    |
| <b>Form/Appearance</b>     | Liquid                                                                         |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                  |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | IgG                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q92947</a>                                                         |
| <b>Dilution Range</b>      | WB:1:500-5000, IHC-P:1:20-200, IP:1:200-2000                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |



Western blot analysis of Mouse kidney(la...



Immunohistochemical staining of human li...



Immunoprecipitation analysis of Hela who...