


Product Datasheet

Human TNFA protein (Active) (orb359192)

Description	Recombinant human TNFA active protein
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His protein
Protein Sequence	MHHHHHH+VRSSSRTPSDKPVAVVNPQAEGQLQWLNRRANALLANGVELR DNQLVVPSEGLYLIYSQVLFKGGCPSTHVLLTHTISRIAVSYQTKVNLLSAIKSPC QRETPEGAEAKPWYEPIYLGGVFQLEKGDRLSAEINRPDYLDFAESGQVYFGIAL
Purity	> 97% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	18.3 kDa
Uniprot ID	P01375
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TNFA
protein ...

