

Product Datasheet

Human TR11B protein (Active) (orb359217)

Description	Recombinant human TR11B active protein
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20 mM PB, pH 6.0, 150 mM NaCl, 0.02 % Tween-80
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Fc-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	ETFPKYLHYDEETSHQLLCDKCPPGYLKQHCTAKWKTVCAPCPDHYYTDSWH TSDECLYCSPVCKELQYVKQECNRTHNRVCECKEGRYLEIEFCLKHRSCPPGFGV VQAGTPERN TVCKRCPDGFSSNETSSKAPCRKHTNCSVFGLLLTQKGNATHDNI CSGNSESTQKCGIDVTL+EPKSSDKTHTCPPCPAPEFEGAPSVFLFPPKPKDTLMI SRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVL TVLHQDWLNGKEYKCKVSNKALPTPIEKTKAKGQPREPQVYTLPPSRDELTKN QVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSR WQQGNVFSCSVMHEALHNHYTQKSLSLSPGK
Purity	> 95% as determined by SDS-PAGE and HPLC.
Source	Yeast
MW	109.6 kDa
Uniprot ID	O00300
Expiration Date	6 months from date of receipt.