



Product Datasheet

H4-K20 antibody (orb35951)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H4C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	11367
Uniprot ID	P62805
NCBI	NP_001029249.1, NP_003536.1, NP_778224.1, NP_003534.1, NP_003529.1, NP_003535.1, NP_003533.1, NP_003530.1, NP_068803.1, NP_003531.1, NP_003532.1, NP_003537.1, NP_003486.1, NP_003539.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse NIH-3T3 c...