

Product Datasheet

Inosine-5'-monophosphate dehydrogenase 2 Polyclonal Antibody (orb359548)

Description	Inosine-5'-monophosphate dehydrogenase 2 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P12268
Expiration Date	12 months from date of receipt.