

Product Datasheet

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 Polyclonal Antibody (orb359592)

Description	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	095167
Expiration Date	12 months from date of receipt.