

Product Datasheet

Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 Antibody (orb359638)

Description	Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y3Z3
Expiration Date	12 months from date of receipt.