

Product Datasheet

CALCR antibody (orb35979)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Synthetic Peptide

Target CALCR

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 55329

Uniprot ID [P30988](#)


NCBI [NP_001158210.1](#), [NP_001158209.1](#), [NP_001733.1](#)

Dilution Range WB: 1:1000, FACS: 1:10-50, IF/ICC: 1:10-50

Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...



Western blot analysis of MCF-7, A431 cel...