




Product Datasheet

AFTPH antibody (orb36062)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	AFTPH
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	102113
Uniprot ID	Q6ULP2
NCBI	NP_982261.2 , NP_001002243.1 , NP_060127.3
Dilution Range	IHC-P: 1:50-100, FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of CEM, Jurkat, MC...