

Product Datasheet

IL1 alpha antibody (orb360764)

Description

IL1 alpha antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human IL1 alpha. Please contact us for the exact immunogen sequence. The peptide is available as orb374974.

Target

IL1 alpha

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

31 kDa

Uniprot ID

[P16598](#), [P01583](#), [Q61730](#)

NCBI

[000576](#), [000567](#)

Entrez


[3553](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:100-2000

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of rat brain...



Immunohistochemical staining of paraffin...



IHC-P image of rat skin tissue using ant...