



Product Datasheet

IL4 antibody (orb360766)


Description	IL4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IL4. Please contact us for the exact immunogen sequence. The peptide is available as orb374976.
Target	IL4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	17 kDa
Uniprot ID	P20096, P07750, P05112
NCBI	000409, 000418
Entrez	3566
Dilution Range	IHC-P: 1:100-500, WB: 1:100-2000
Expiration Date	12 months from date of receipt.



IHC-P image of human thyroid tumour tissue...



IHC-P staining of rat thymus tissue using...



Western blot analysis of rat lymph node

...