

Product Datasheet

Myeloperoxidase antibody (orb360767)

Description

Myeloperoxidase antibody


Species/Host


Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated


Tested

IHC-P

Applications

Immunohistochemical staining of rat sple...

Immunogen

KLH conjugated synthetic peptide derived from human Myeloperoxidase. Please contact us for the exact immunogen sequence. The peptide is available as orb374977.

Target


Myeloperoxidase

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml



Immunohistochemical staining of paraffin...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

84 kDa

Uniprot ID

P11247, **P05164**, **D3ZGE2**

NCBI

000250, **000241**

Entrez

4353

Dilution Range

IHC-P: 1:100-500

Expiration Date

12 months from date of receipt.