



Product Datasheet

CLU antibody (orb360771)


Description	CLU antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic peptide (please contact us for details)
Target	CLU
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Monoclonal antibodies are purified by peptide affinity chromatography
MW	52 kDa
Uniprot ID	P05371, P10909, Q06890
NCBI	001822, 001831
Entrez	1191
Dilution Range	IHC-P: 1:100-500, WB: 1:100-3000
Expiration Date	12 months from date of receipt.



WB analysis of Clu recombinant protein 1...



IHC-P image of human breast cancer tissu...



WB analysis of Clu recombinant protein 2...