

Product Datasheet

Apolipoprotein CIII antibody (FITC) (orb360850)

Description	Apolipoprotein CIII antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Recombinant protein
Target	Apolipoprotein CIII (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Uniprot ID	P02656 , P33622 , P06759
NCBI	NM_000040.1 , NP_000031.1
Entrez	345
Dilution Range	IF/ICC: 1:30-300
Expiration Date	12 months from date of receipt.