

Product Datasheet

ACSF2 antibody (FITC) (orb36087)

Description	ACSF2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Acyl-CoA synthetase family member 2, mitochondrial protein (42-615AA)
Target	ACSF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96CM8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.