

Product Datasheet

TSPYL4 antibody (orb36116)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TSPYL4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45126
Uniprot ID	Q9UJ04
NCBI	NP_067680.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

