


Product Datasheet

SLC7A3 antibody (orb36232)

Description	SLC7A3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Cationic amino acid transporter 3 protein (428-475AA)
Target	SLC7A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WY07
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of MCF-7 cell...