

Product Datasheet

SARS virus Sn antibody (orb36283)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	S
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	139125
Uniprot ID	P59594
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

