

Product Datasheet

METTL2A antibody (FITC) (orb36292)

Description	METTL2A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Methyltransferase-like protein 2A protein (1-242AA)
Target	METTL2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96IZ6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.