



Product Datasheet

HAND2 antibody (orb36304)

Description	Rabbit polyclonal antibody to HAND2 (also named as Thing-2 and Hed, dHAND) is a 24 kDa, class B memb...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	HAND2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23666
Uniprot ID	P61296
NCBI	NP_068808.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of HAND2 antibody ...