




Product Datasheet

ICOS antibody (orb36326)

Description	Rabbit polyclonal antibody to ICOS is highly expressed on tonsillar T-cells. It promotes all basic T...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ICOS
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22625
Uniprot ID	Q9Y6W8
NCBI	NP_036224.1
Dilution Range	IHC-P: 1:50-100, WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of Jurkat cells...