




Product Datasheet

Leptin Receptor antibody (orb36330)

Description	Rabbit polyclonal antibody to Leptin Receptor
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	LEPR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	132494
Uniprot ID	P48357
NCBI	NP_001003679.1 , NP_001185618.1 , NP_001185617.1 , NP_002294.2 , NP_001003680.1 , NP_001185616.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of K562 cells (...)