

Product Datasheet

LY 344864 (orb363706)

Description	LY344864 is a selective receptor agonist with an affinity of 6 nM (Ki) at the recently cloned 5-HT1F...
Note	For research use only
Formula	C ₂₁ H ₂₂ FN ₃ O
MW	351.4173
CAS Number	[186544-26-3]
SMILES	<chem>O=C(C1=CC=C(F)C=C1)NC2=CC=C3NC(CC4)=C(C[C@@H]4N(C)C)C3=C2</chem>
Expiration Date	12 months from date of receipt.