




Product Datasheet

NRG1 antibody (orb36387)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	Synthetic Peptide
Target	NRG1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70392
Uniprot ID	Q02297
NCBI	NP_039254.1 , NP_039256.2 , NP_001153471.1 , NP_039258.1 , NP_039250.2 , NP_039252.2 , NP_001153480.1 , NP_001153467.1 , NP_039253.1 , NP_001153477.1 , NP_039251.2 , NP_001153473.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
NRG1
antibod...



Western blot
analysis of
anti-NRG1
antib...