

Product Datasheet

(R)-1,2,3,4-Tetrahydro-3-isoquinolinecarboxylic acid (orb363924)

Description	(R)-1,2,3,4-Tetrahydro-3-isoquinolinecarboxylic acid is a constrained Phe analogue which can fold in...
Note	For research use only
Formula	C ₁₀ H ₁₁ NO ₂
Purity	> 98%,Standard References
MW	177.1998
CAS Number	[103733-65-9]
SMILES	<chem>O=C([C@@H]1NCC2=C(C=CC=C2)C1)O</chem>
Expiration Date	12 months from date of receipt.