

Product Datasheet

SR9009 (orb363935)

Description	SR9009 is an agonist of REV-ERB with IC ₅₀ = 670 nM for REV-ERB α and IC ₅₀ = 800 nM for REV-ERB β .
Note	For research use only
Formula	C ₂₀ H ₂₄ ClN ₃ O ₄ S
MW	437.9402
CAS Number	[1379686-30-2]
SMILES	CCOC(=O)N1CCC(C1)CN(CC2=CC=C(C=C2)Cl)CC3=CC=C(S3)[N+](=O)[O-]
Expiration Date	12 months from date of receipt.