

Product Datasheet

PI-1840 (orb364002)

Description	PI-1840 is a potent and selective inhibitor for chymotrypsin-like (CT-L) (IC50 value = 27 ± 0.14 nM...
Note	For research use only
Formula	$C_{22}H_{26}N_4O_3$
MW	394.4668
CAS Number	[1401223-22-0]
SMILES	<chem>O=C(N(C(C)C)CC1=NC(C2=CC=CN=C2)=NO1)COC3=CC=C(CCC)C=C3</chem>
Expiration Date	12 months from date of receipt.