

Product Datasheet

Ingenol (orb364036)

Description	Ingenol is an extremely weak protein kinase C (PKC) activator.
Note	For research use only
Formula	C ₂₀ H ₂₈ O ₅
MW	348.4333
CAS Number	[30220-46-3]
SMILES	<chem>O([H])C12C([H])(C(C([H])([H])[H])=C([H])C31C(C([H])(C([H])=C(C([H])([H])O[H])C2([H])O[H])C1([H])C(C([H])([H])[H])(C([H])([H])([H])C1([H])C([H])([H])C3([H])C([H])([H])[H])=O)O[H])</chem>
Expiration Date	12 months from date of receipt.