

Product Datasheet

[APLILSR]-pPSA (orb364054)

Description	Peptide to [APLILSR]-pPSA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Pro-Leu-Ile-Leu-Ser-Arg-Ile-Val-Gly-Gly-Trp-Glu-Cys-Glu-Lys
Formula	$C_{80}H_{131}N_{21}O_{22}S$
Purity	> 95%
Source	Synthetic
MW	1771.12
Expiration Date	12 months from date of receipt.