

Product Datasheet

[Asn370]-Tyrosinase (368-376) (orb364060)

Description	Peptide to [Asn370]-Tyrosinase (368-376)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Met-Asn-Gly-Thr-Met-Ser-Gln-Val
Formula	$C_{42}H_{67}N_{11}O_{15}S_2$
Purity	> 95%
Source	Synthetic
MW	1030.19
Expiration Date	12 months from date of receipt.