

Product Datasheet

[Des-Arg30, Des-Pro31]-Dendroaspis Natriuretic Peptide (orb364069)

Description	Peptide to [Des-Arg30, Des-Pro31]-Dendroaspis Natriuretic Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Val-Lys-Tyr-Asp-Pro-Cys-Phe-Gly-His-Lys-Ile-Asp-Arg-Ile-Asn-His-Val-Ser-Asn-Leu-Gly-Cys-Pro-Ser-Leu-Arg-Asp-Pro-Asn-Ala-Pro-Ser-Thr-Ser-Ala (Disulfide bridge Cys7-Cys23)
Formula	$C_{169}H_{263}N_{51}O_{54}S_2$
Purity	> 95%
Source	Synthetic
MW	3937.42
Expiration Date	12 months from date of receipt.