

Product Datasheet

[Gln11]-Beta-Amyloid (1-40) (orb364078)

Description	Peptide to [Gln11]-Beta-Amyloid (1-40)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Gln-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val
Formula	$C_{194}H_{296}N_{54}O_{57}S$
Purity	> 95%
Source	Synthetic
MW	4328.91
Expiration Date	12 months from date of receipt.