



Product Datasheet

HSC70 Interacting Protein HIP antibody (orb36408)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ST13
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41332
Uniprot ID	P50502
NCBI	NP_003923.2
Dilution Range	IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in A2058 cell l...