

Product Datasheet

[Gln22]-Beta-Amyloid (6-40) (orb364080)

Description	Peptide to [Gln22]-Beta-Amyloid (6-40)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Gln-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val
Formula	$C_{167}H_{258}N_{46}O_{48}S$
Purity	> 95%
Source	Synthetic
MW	3710.26
Expiration Date	12 months from date of receipt.