

Product Datasheet

[His11]-Substance P (orb364090)

Description	Peptide to [His11]-Substance P
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-His-NH ₂
Formula	C ₆₄ H ₉₆ N ₂₀ O ₁₃
Purity	> 95%
Source	Synthetic
MW	1353.60
Expiration Date	12 months from date of receipt.