

Product Datasheet

[pY185]-MAP (177-189) kinase (orb364106)

Description	Peptide to [pY185]-MAP (177-189) kinase
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-His-Thr-Gly-Phe-Leu-Thr-Glu-pTyr-Val-Ala-Thr-Arg
Formula	$C_{67}H_{100}N_{18}O_{22}$
Purity	> 95%
Source	Synthetic
MW	1509.65
Expiration Date	12 months from date of receipt.