

Product Datasheet

[Tyr1]-TRAP-7 (orb364119)

Description	Peptide to [Tyr1]-TRAP-7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Phe-Leu-Leu-Arg-Asn-Pro
Formula	$C_{455}H_{67}N_{11}O_{10}$
Purity	> 95%
Source	Synthetic
MW	922.10
Expiration Date	12 months from date of receipt.