

Product Datasheet

234 CM (orb364122)

Description	Peptide to 234 CM
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Tyr-Ile-Cys-Asn-Ser-Ser-Cys-Met
Formula	$C_{42}H_{69}N_{11}O_{14}S_3$
Purity	> 95%
Source	Synthetic
MW	1048.26
Expiration Date	12 months from date of receipt.