

Product Datasheet

2A/2B Dengue Protease Substrate (orb364127)

Description	Peptide to 2A/2B Dengue Protease Substrate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Arg-Thr-Ser-Lys-Lys-Arg- pNA
Formula	$C_{39}H_{68}N_{16}O_{11}$
Purity	> 95%
Source	Synthetic
MW	937.08
Expiration Date	12 months from date of receipt.