

Product Datasheet

2B/3, Dengue Protease Substrate (orb364129)

Description	Peptide to 2B/3, Dengue Protease Substrate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Glu-Val-Lys-Lys-Gln-Arg- pNA
Formula	$C_{41}H_{58}N_{14}O_{12}$
Purity	> 95%
Source	Synthetic
MW	949.09
Expiration Date	12 months from date of receipt.